ADHD Catcher

Clinician Operations Manual · 임상가 운영 매뉴얼

 Version 11.0 — Adaptive Treatment Engine

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ENGLISH

ADHD Catcher v11.0 — Adaptive Treatment Engine is a multi-engine, learning-based decision support system for neurofeedback / learning-performance clinics.

It does NOT:

- Make medical diagnoses
- Replace clinical judgment
- Prescribe medications
- Make treatment decisions autonomously

It DOES:

- Analyze QEEG/HRV/survey data with normative comparison
- Track patient progress over sessions (trajectory analysis)
- Compare with synthetic + real cohorts (continuously learning)
- Recommend next-step decisions with confidence scores
- Provide explainable AI reasoning (XAI) for every recommendation
- Detect critical safety factors automatically
- Document every decision for medical-legal protection
- Support clinicians — final authority remains with clinician

v11.0 Key Innovations:

-  **5 Core Math Functions** — formal mathematical foundation
-  **3-Stage System** — Engagement → Validation → Integration
-  **Threshold-Free Architecture** — all numbers in policies, not code
-  **Multi-Tag Policy System** — flexible policy combination
-  **Self-Correcting Synthetic DB** — improves with every patient
-  **Critical Factor Detection** — automatic safety monitoring

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v11

 ENGLISH

v11.0 represents a paradigm shift from "static analysis tool" to "continuously learning treatment engine."

Core Principles

- The system has a brain.** Not just a calculator. The 5 core math functions form a coherent decision-making system inspired by AlphaGo's evaluation paradigm — but adapted for the 3D space of clinical neurofeedback.
- 3D Decision Space.** Every decision integrates: (1) State (current patient vector), (2) Time (trajectory & velocity), (3) Context (BPS + goals + capacity).
- Progressive measurement.** Different metrics matter at different stages: subjective experience first (engagement), objective data middle (validation), functional outcome last (integration).
- Hardcoded numbers are dead code.** All thresholds, weights, and rules live in editable policies. The engine itself never changes — policies evolve.
- AI assists, clinician decides.** Every recommendation includes confidence, reasoning, and alternatives. Clinician can always override (and overrides become learning data).
- Critical factors get specialized treatment.** Suicide risk, dangerous substance combinations (e.g., benzo + alcohol), active psychosis — these never go through the standard pattern-matching algorithm. They route to specialized safety protocols.
- Self-improving system.** Synthetic cohort starts the cold-start problem. Real patient data progressively replaces synthetic. The system gets smarter every session.

What Changed from v3.0

Aspect	v3.0 (Old)	v11.0 (New)
Decision basis	Rule-based	Math-function based + policy
Time tracking	None	Full trajectory analysis
Cohort comparison	Static reference	Self-evolving synthetic + real
Stage awareness	None	3-stage system with transitions
Policy	Hardcoded thresholds	Multi-tag editable policies
Safety	Manual review	Auto critical factor detection
Documentation	Manual notes	Auto audit trail of every decision
Learning	None	Continuous (every patient)

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3. 5 Core Math Functions · 5개 핵심 수학 함수 ★ NEW v11

EN: The mathematical heart of v11.0. Inspired by AlphaGo's evaluation function but adapted for 3D clinical decision space (State × Time × Context).

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$$\text{Decision} = \pi(S, T, C)$$

EN: Where S = State, T = Time, C = Context. The policy function π returns: action + confidence + reasoning.

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1 S(patient, t) — State Function · 상태 함수

🇺🇸 ENGLISH

Represents patient as a 13-dimensional vector at time t:

- **QEEG (5):** theta/beta, frontal alpha asymmetry, alpha peak, posterior alpha, coherence
- **Scales (3):** primary, secondary, BPS total
- **Function (2):** goal achievement %, daily function
- **Subjective (1):** patient self-report
- **Observation (1):** clinician observation
- **Time (1):** session number

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2 V(S, C) — Value Function · 가치 함수

EN: "How good is this state for this patient?" Returns 0 (worst) to 1 (perfect). Equivalent to AlphaGo's board evaluation.

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EN: Weights change by stage:

- 🌱 Initial: subjective + observation dominate
- 🌱 Active: QEEG + scales dominate
- 🌱 Integration: goal achievement dominates

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3 τ (patient, lookback) — Trajectory Function · 궤적 함수

EN: Analyzes how patient is changing over recent N sessions. Returns:

- **Direction:** Where is patient heading?

- **Velocity:** How fast?
- **Acceleration:** Speeding up or slowing down?
- **Volatility:** Stable or oscillating?
- **Plateau detection:** Has progress stalled?
- **Confidence:** How reliable is this trajectory?

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4 Q(S, A, C) — Action Value Function · 행동 가치 함수

EN: "If we do action A in state S with context C, what's the expected outcome?" Learned from synthetic cohort + real patient outcomes.

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EN: Returns:

- **Expected value:** 0-1 score
- **Confidence:** Based on similar case count
- **Sample size:** How many comparable cases
- **Distribution:** Range of outcomes (for uncertainty visualization)

5 $\pi(S, C)$ — Policy Function · 정책 함수 ★

EN: The final decision function. Combines all of the above. This is what gets called every session.

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EN — Step-by-step:

1. Evaluate Q(S, A, C) for all possible actions
2. Get trajectory τ for context
3. Filter actions through critical factor safety check
4. Select best action with exploration-exploitation balance
5. Generate explainable reasoning
6. Predict expected future trajectory
7. Quantify uncertainty (epistemic + aleatoric)

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4. Stage System · 단계 시스템 ★ NEW v11

EN: Treatment progresses through 3 stages with different evaluation priorities and clinical goals. Stage transitions are hybrid (smart auto + safety guardrails + clinician override).

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Stage 1: INITIAL (Engagement Phase)

EN — Sessions 1-8 (typical): Build patient trust and establish baseline. The patient should feel something is happening — even if measurements haven't moved yet.

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Stage 2: ACTIVE (Validation Phase)

EN — Sessions 9-20 (typical): Establish objective evidence of change. QEEG and validated scales should show measurable improvement.

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Stage 3: INTEGRATION (Outcome Phase)

EN — Sessions 21+ (typical): Generalize gains to daily life. Functional outcomes (school, work, relationships) become primary measure.

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Stage Transition Logic · 단계 전환 로직

ENGLISH

Hybrid system — three pathways to advance:

-  **Early Advance (data-driven):** Patient is responding faster than expected. System auto-suggests advancement.
-  **Standard Advance (session-based):** Standard expected progression timeline.
-  **Late Advance (safety guardrail):** Patient is slower but maximum sessions reached. Forced advance with monitoring.

Special triggers:

-  **Plateau detection:** Auto-advance to integration if QEEG changes < 0.2 SD over last 5 sessions (despite incomplete advance criteria)

- 🚨 **Critical factor:** Override standard transitions, route to specialized algorithm
- 🎓 **Graduation:** Goal achievement $\geq 85\%$ + 4+ maintenance sessions → recommend course completion

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URL

Production URL:

https://neurocatchers.com/neurocatchers_multifile/catchers/adhd-catcher.html

Admin Rules Editor: <https://neurocatchers.com/admin/>

Browser Requirements · 브라우저 요구사항

🇺🇸 ENGLISH

- Chrome / Edge / Safari / Firefox (latest 2 versions)
- JavaScript enabled
- Minimum 1280x800 screen resolution recommended
- localStorage enabled (stores patient data locally)

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⚠️ Important · 중요:

EN: Data is stored in browser localStorage. Use a dedicated clinic computer. Export data regularly as backup.

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~1 min

Click "+ New" in Patient panel (top left).

Enter: Name, Age, Gender, any relevant notes.

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~2 min

Goal (0-100): Bio / Psycho / Social targets based on clinical judgment.

Protocol: Choose from SMR / Beta / Theta_inhibit / Alpha_theta / ILF / HRV_BF.

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~30-45 min

Record QEEG (19-channel, Eyes Open 5min + Eyes Closed 5min + cognitive task) → upload via Mitsar.

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Manual Input: Enter key values directly.

Upload File: JSON / CSV / Mitsar / Neuroguide export.

Watch *Signal Inventory* go green — aim for 80%+ completeness.

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<10 sec

Click " Run Brain Engine".

Layers 0-5 auto-analyze: Topographic → Features → States → 3-Axis → Profile → Clinical Impression.

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~10 sec

Click " Save Session" in Layer 3.5 (3D Spatial Trajectory).

Session is saved to patient history for trajectory tracking.

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7 Review & discuss · 검토 및 상담

15-20 min

Show patient: Topographic → 3-Axis → Clinical Impression.

Explain in plain language. Schedule next session.

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Trajectory, Cohort, Decision, Weight 확인.

4 Clinical decision · 임상 결정

Use Decision Engine recommendation as **input**, not command.

Integrate with clinical judgment. Adjust weights if patient/context requires.

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7. UI Walkthrough · 화면 설명 (with Decision Cockpit)

Left Panel · 왼쪽 패널

Section	Function · 기능
👤 Patient	Create/switch patients · 환자 생성/전환, Export/Import JSON

Section	Function · 기능
 Signal Inventory	Data completeness tracker · 데이터 완성도 (EEG/HRV/Survey/Cognitive/Context)
 Actions	Manual Input / Surveys / Upload File / Sample cases
 Run Brain Engine	Orange/green button — runs full analysis · 전체 분석 실행
 Session Timeline	Per-patient session history · 환자별 세션 히스토리
 Goal & Protocol	Set target BPS + protocol selection · 목표 BPS 및 프로토콜 선택

Right Panel (Layers) · 오른쪽 패널 (레이어)

Layer	Content · 내용	Activated by
0  Brain Topographic	19-channel QEEG heatmap · 뇌파 공간 분포	Run Brain Engine
1 Feature Extraction	Z-scores vs normative DB · Normative DB와 비교	Run Brain Engine
2 State Inference	Cortical arousal, attention, etc. · 각성도, 주의 등	Run Brain Engine
3 3-Axis Mapping	Bio / Psycho / Social scores	Run Brain Engine
3.5  3D Spatial Trajectory	3D position + Save Session	Run Brain Engine
4 Profile Matching	Inattentive/Hyperactive/Combined %	Run Brain Engine
5 Clinical Impression	Narrative + recommendations	Run Brain Engine
 6 Trajectory Engine	Velocity/Alignment/ETA + 4 charts	Need ≥2 sessions
 7 Cohort Comparison	500 synthetic patients · Percentile	Click "Generate Cohort"
 8 Adaptive Decision	Auto-generated recommendation	Automatic after Layer 6
 9 Weight Adjustments	Modal with sliders for 14 weights	Click "Adjust Weights"

🎯 Decision Cockpit · 결정 콕핏 ★ NEW v11

EN: The Decision Cockpit appears after Layer 8 and provides a high-level summary of the AI's recommendation. This is the primary interface for clinical decision-making.

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Example Cockpit View · 콕핏 화면 예시

 Recommended Action: **MODIFY_Adjust**
Confidence: 78% 

 Why?

- Current state value: 0.42 (mid-range)
- Recent trajectory: improving (+0.05/session)

- 23 similar patients improved with this action
- Expected outcome: 0.65 (range 0.55-0.75)

 Alternatives considered:

- MAINTAIN: 0.48 (-0.17 vs MODIFY) 
- SWITCH: 0.62 (-0.03, similar) 
- PAUSE: 0.30 (not recommended) 

 Uncertainties:

- Patient has comorbid anxiety (data weaker)
- Only 23 similar cases (low confidence)

 Evidence: Arns 2014 + BNM cohort 87% match

[ Accept] [ Modify] [ See full analysis]

Override Workflow · Override 워크플로

EN: When clinician overrides AI recommendation:

1. System asks for override reason (free text or category)
2. Decision is recorded in patient record + audit log
3. Override pattern is tracked for system learning
4. Clinician's chosen action is implemented

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EN: Decision Engine uses trajectory + cohort + critical factors to recommend one of 9 decisions. Use as input, not command.

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🏥 REFER_Medical

EN: Medical concern — refer to physician.
KR: 의학적 우려 — 의사 의뢰.

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Decision Priority Order · 결정 우선순위

1. **Medical concern** (critical) → REFER_Medical
2. **Goal reached** (distance < 3.0) → STOP
3. **Reversing** (alignment < -0.3) → MODIFY_Switch
4. **Alternative diagnosis** → REFER_Diagnostic
5. **Plateau + confounders** (sleep/stress/med) → PAUSE
6. **Cohort bottom 25%** (after ≥5 sessions) → MODIFY_Switch
7. **Plateau** (after ≥5 sessions) → MODIFY_Adjust
8. **Drifting** (after ≥5 sessions) → MODIFY_Add
9. **Excellent progress** → MAINTAIN (high confidence)
10. **Default (moderate progress)** → MAINTAIN (medium confidence)

9. Policy System · 정책 시스템 ★ NEW v11

EN: All thresholds, weights, and rules live in **editable policies** — not in code. This is the "Threshold-Free Architecture" principle.

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EN: Each policy has tags. The system auto-combines policies based on patient profile.

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Category · 분류	Examples · 예시
 Goal-based · 목표 기반	focus_enhancement, anxiety_reduction, sleep_optimization, peak_performance ★ (BNM specialty), learning_enhancement ★ (BNM specialty)
 Diagnosis-based · 진단 기반	adhd_inattentive, adhd_combined, gad, mdd, ptsd
 Demographic · 인구학적	child_6_12, teen_13_17, adult_18_64, senior_65_plus
 Style · 치료 스타일	aggressive (fast), standard, conservative (safety-first)
 Overlay · 덧붙임	anxiety_comorbid, medication_active, trauma_informed, substance_history

Example Policy Selection · 정책 선택 예시

EN — Patient: 32F, ADHD-Inattentive + GAD + cognitive enhancement goal
System auto-suggests:

- Adult Attention Enhancement (87% match)
- + Anxiety Co-treatment Overlay
- + Female 30s Cognitive Performance overlay

Combined policy preview:

- QEEG threshold: 1.5×1.2 (anxiety modifier) = 1.8
- Scales monitored: ASRS + GAD-7
- Stage transitions: standard + HRV stability check
- Estimated course: 24 sessions

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 **Conflict Resolution · 충돌 해결**

EN: When policies conflict (e.g., aggressive style says "advance fast" but anxiety overlay says "go slow"), the system uses **Hybrid resolution**:

- **Default:** Choose the safer/more conservative value (medical safety first)
- **Override:** Clinician can explicitly choose differently
- **Logging:** All overrides are recorded for system learning

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v11

EN: Solves the cold-start problem of real-world NF data scarcity by generating literature-based virtual patients. Self-corrects as real patient data accumulates.

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Hybrid Approach (Literature-Anchored Simulation)

EN:

1. **Literature anchors** — Mean \pm SD from meta-analyses (e.g., Arns 2014 n=600)
2. **Patient simulation** — Generate 500 virtual patients with realistic distributions
3. **Trajectory generation** — Each patient simulated session-by-session with noise + drift
4. **Responder categories** — Fast (20%), Standard (45%), Slow (20%), Non-responder (15%)
5. **Confounder application** — Age, comorbid anxiety, medication, adherence affect outcomes

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EN: As real patient data flows in:

Real Patients	Synthetic Weight	Real Weight
0-30	95%	5%
30-100	60%	40%
100-300	20%	80%
300+	5%	95% (synthetic graduated)

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Validation Score · 검증 점수

EN: System dashboard shows overall validity (0-100):

- Mean match (20%) · Distribution match (25%) · Correlation match (20%) · Trajectory match (20%) · Responder rate (15%)
- **EXCELLENT** ≥ 90 | **GOOD** ≥ 75 | **ACCEPTABLE** ≥ 60 | **NEEDS REVIEW** < 60

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EN: Certain conditions require **specialized algorithms**, NOT general pattern matching. The system auto-detects and routes appropriately.

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Tier 3 — Specialized Algorithm Required

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Tier 3 Alert Workflow · Tier 3 경고 워크플로

EN: When Tier 3 detected, the system:

1.  Shows critical alert modal
2.  Lists detected factors + risks
3.  Provides system recommendations (defer/refer/specific actions)
4.  Requires clinician decision: Accept / Override / Custom plan
5.  If override, requires reason (recorded for audit)
6.  Auto-generates documentation (medical-legal protection)

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12. Weight Adjustments · 가중치 조정

ENGLISH

Adjust how the Decision Engine weighs different factors. Useful when:

- Patient has specific history (e.g., trauma)
- Patient preferences (e.g., avoid stimulants)
- Unusual presentation needs custom weighting

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Tip · 팁:

EN: Every adjustment is logged with timestamp + reason. Document your reasoning in the adjustment reason field for compliance.

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All data is in browser `localStorage` :

- `nc_patients` — Patient records · 환자 기록
- `nc_sessions_adhd` — ADHD sessions · ADHD 세션
- `nc_active_patient_adhd` — Active patient · 활성 환자

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⚠ Critical · 중요:

EN: Clearing browser data = PERMANENT data loss. Always backup before browser maintenance.

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15. Clinical Safety Notes · 임상 안전 주의사항

REFER_Medical triggers · REFER_Medical 발동 시

EN: If system flags medical concern, stop NF training and refer to physician. Document all findings.

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Non-physician Scope · 비의사 범위

ENGLISH

- Use "training/coaching/performance optimization" language
- Avoid "treatment/cure/diagnosis" terminology
- Refer out for medical/psychiatric diagnosis
- Document clinical judgment and referral rationale
- Complies with MA non-physician practitioner scope

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EN: Every session should include: data sources, analysis summary, decision rationale, any weight adjustments with reasons, next steps.

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ADHD Catcher v3.0 · Clinician Manual

Boston Neuromind LLC · Canton, Massachusetts

© 2026 · Confidential · Clinic internal use only

This tool supports clinical decision-making but does not replace clinical judgment.

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