# ADHD Catcher

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ADHD Catcher가 뭐예요?

## ENGLISH

**ADHD Catcher** is a smart brain-analysis tool used at our clinic.

It takes your QEEG (brain wave), HRV (heart rhythm), and survey data to:

- ✨ See how your brain is working right now
- ✨ Compare with healthy patterns
- ✨ Track how you're improving over time
- ✨ Help your clinician choose the best training plan for you

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## Tip · 팁

**EN:** You don't need to understand the technical details! The clinician will explain everything.

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### For patients — Very simple! (4 steps)

#### 1 Relax — we do the work

You just sit comfortably while we record your brain waves (QEEG) and heart rhythm (HRV).

*Takes about 20-30 minutes.*

#### 2 Answer a few questions

Fill in some simple surveys about how you've been feeling — sleep, focus, mood.

*10-15 minutes.*

#### 3 Clinician runs the analysis

Our clinician clicks "**Run Brain Engine**" — the system analyzes everything in seconds.

You get a colorful report showing what's going on.

#### 4 Discuss & plan together

The clinician walks you through the results in plain language.

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### 3. Your Treatment Journey

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#### ENGLISH

Your treatment progresses through **3 stages**. At each stage, we focus on different signs of progress so you can see real change at every point.

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#### ENGLISH

##### Sessions 1-8 (typical)

- 🌟  Goal: **You feel something is changing**
- 🌟  We focus on: How **you** feel and what **we observe** in session
- 🌟  Why: Trust + early signals build motivation to continue
- 🌟  Sessions can be done at your own pace (1-3x per week)

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#### ENGLISH

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### Sessions 9-20 (typical)

🌟 🎯 Goal: **Objective evidence of brain change**

🌟 📊 We focus on: **QEEG measurements + standard scale scores**

🌟 💡 Why: This is when scientific change becomes measurable

🌟 📈 You'll see clearer numbers and graphs

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### 🇺🇸 ENGLISH

#### Sessions 21+ (typical)

🌟 🎯 Goal: **Real-life improvement** (school, work, relationships)

🌟 📊 We focus on: **Goal achievement** in your daily life

🌟 💡 Why: This is the real reason you came — to improve life

🌟 🎓 When goals are reached, we recommend graduation

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### 💡 Important · 중요

**EN:** Sessions are not "one-time use." We track your progress across all sessions to give you smarter recommendations as we learn more about your brain. The system gets better at helping you with every session.

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## 4. What you will see · 무엇이 보이나요?

**EN:** You'll see 10 colorful panels. Don't worry — you don't need to understand all of them. Here are the important ones.

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#### ENGLISH

**What you see:** A colorful "heatmap" of your head showing which brain regions are more or less active.

**What it means:** Red = high activity, Blue = low activity. Your clinician will interpret this!

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### 3-Axis Mapping · 3축 지도 (Bio · Psycho · Social)

#### ENGLISH

Three numbers (0-100) showing three aspects of your wellness:

- ✨ **Bio** — Body (sleep, HRV, arousal)
- ✨ **Psycho** — Mind (attention, focus, emotion)
- ✨ **Social** — Daily life (relationships, school/work)

Higher = better. Goal is to move these up over time.

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### ENGLISH

After 2+ sessions, you'll see:

-   Charts showing your improvement
-   How close you are to your goal
-   Estimated sessions until goal
-   Whether you're on track

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### ENGLISH

The system compares you with **500 similar patients** (computer-generated based on research).

You'll see: "Top 25%" or "Average" or "Below average" — helping you understand your progress.

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### Q. How often should I come? · 얼마나 자주 와야 해요?

**EN:** Usually every 2-4 weeks, depending on your plan. Your clinician will advise.

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**EN:** Yes! All data is stored securely in our clinic system. You can request export at any time.

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### Q. Do I need a diagnosis to use this? · 진단이 있어야 해요?

**EN:** No. This tool is for performance/wellness training. It's not a medical diagnosis tool.

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**EN:** Usually 5-15 sessions, depending on protocol. The Trajectory chart helps you track!

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**EN:** Yes! The system adapts to children, teens, and adults. Age-specific comparisons are used.

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**EN:** No problem! Your progress is saved. Just resume when you come back.

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## 7. Questions? · 문의

 **ENGLISH**

**Your clinician is your first contact!**

For any questions about your results, treatment plan, or the system — please ask during your session.

**Boston Neuromind LLC**  
Canton, Massachusetts

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**Boston Neuromind LLC**  
매사추세츠주 캔톤

*This is a performance/wellness training tool, not a medical diagnosis system.*

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